

Resource Material – Speech and Language

Category	Type	Title	Annotation
MF000013	<i>Adaptive/Assistive Technology</i>	<i>Super Hawk</i>	The Superhawk has a highly sensitive touchpad. It also has alterable dimensions, level changing, message buffer, and game commands. Superhawk also has auditory scanning with options: digital speck, 4 single switch inputs and 4 relay outputs. New features include an LCD screen, a view mode to see what is programmed into a page or cell and see how much memory is available. It also has upgradeable memory.
MF 00014	<i>Adaptive/Assistive Technology</i>	<i>Hawk</i>	The Hawk is an easy and quickly programmable augmentative communications device employing digital voice technology for high quality, natural speaking voice output. A built-in microphone allows the recording of up to 5 seconds worth of digital speech into each of the nine selection areas. The Hawk's weight and hand portability make it suitable for field trip activities.
MF000015	<i>Kit</i>	<i>ASHA Seminar School Reentry Following Traumatic Brain Injury</i>	This document is designed to guide school staff in their development of IEPs for students with TBI.
MF000026	<i>Adaptive/Assistive Technology</i>	<i>Whisper Wolf</i>	Whisper Wolf is a fully operational user programmable MEGA WOLF which also provided auditory scanning and ingle switch selection. It uses synthesized (robotic) voice. For users able to point with approximately 6 ounces of pressure. Includes 18PP. demonstration on vocabulary, complete with picture overlays, customized vocabulary for school, and 2 additional memory areas for programming. Overlays and power adapter.
MF000026	<i>Book</i>	<i>Pocket guide to Assessment in Speech Language Pathology</i>	This is a pocket guide resource book that summarizes basic information on disorders of communication and specifies steps involved in assessing them.
MF000027	<i>Book</i>	<i>Pocket Guide to Treatment in Speech Language Pathology</i>	The pocket guide gives the essences of treatment in a step by step form that promotes easy understanding and ready reference just before beginning treatment.
MF000034	<i>Book</i>	<i>The Linda Mood Phoneme Sequencing Program for Reading, Spelling and Speech (LIPS)</i>	This major revision of the popular Auditory Discrimination in depth (ADD) program gives you the tools to develop the phoneme sequencing ability that research documents as the best single predictor of success in learning to read and spell. The LIPS clinical version includes a trainer's manual, research booklet, overview video tape, audio cassette, 15 colored felts, I set of nose and ear felts, 1 set of 24 colored plastic blocks, 1 set of letter tiles and 1 set of 5 static cling letter-symbols and mouth pictures.